



INJEKT 10ml Next Generation – DANE TECHNICZNE

UWAGA : między oboma modelami strzykawk INJEKT® 10 ml poprzedniej i nowej nie ma różnicy w :

- przeznaczeniu
- funkcjonalności
- obudzie
- materiale
- nazwie produktu

Poprzednia wersja

NOWA wersja

Widok		
Nazwa	INJEKT® 10ml	INJEKT® 10ml
Nr kat.	4606108V	4606108 N
Typ	Strzykawka 2-częściowa	Strzykawka 2-częściowa
Końcówka	Luer niecentryczna	Luer niecentryczna
Skala co	0,5 ml	0,5 ml
Całkowita skala	10 ml +2ml	10 ml (brak rozszerzonej skali)

Komora strzykawki

długość 86,00 mm



10 → 12ml

krótsza o ok 7mm



10 ml

INJEKT 10ml Next Generation – DANE TECHNICZNE

<p>Tłok strzykawki</p>	<p>długość 88,70 mm dłuższy, stabilizujące go w komorze 4 skrzydełka</p>  <p>A photograph of a teal-colored syringe plunger. It has a long, cylindrical body with four distinct wings at the base. A purple label with the word "current" is positioned to the left of the plunger.</p>	<p>krótszy o ok. 8mm "wycieniowany" w kierunku uchwytu, stabilizujące go w komorze 3 skrzydełka</p>  <p>A photograph of a teal-colored syringe plunger. It is shorter than the current model and has a tapered, hourglass-like shape with three wings at the base. A green label with the word "new" is positioned to the right of the plunger.</p>
-------------------------------	--	---